

# See More. Do More.

*With Personal Ultrasound*

**signos**  
Personal Ultrasound

Deliver patient care quickly and affordably with Signos Personal Ultrasound — the world's first palm-sized, high-performance ultrasound system.






### SIGNOS PERSONAL ULTRASOUND—THERE WHEN YOU NEED IT

Weighing less than 300 grams and featuring superb image resolution, Signos is the smallest and most cost-effective ultrasound system available today. Its long battery life and quick start-up time make it always ready for immediate scanning. Signos Personal Ultrasound delivers the best of miniaturization technology right to your practice. **Wear it like a stethoscope. Put it in your pocket. Take it anywhere you need it. Signos is ready when you are.**

### PROCEDURAL AND DIAGNOSTIC CONFIDENCE

Signos delivers critical diagnostic point-of-care information to you—right in the palm of your hand. With a simple, easy-to-use interface and minimal capital investment, Signos delivers profound imaging solutions for point-of-care medicine and practical exams. Operating at either 3.5 MHz or 7.5MHz, Signos facilitates a wide range of clinical applications, including ultrasound guidance for needle placement, AAA exams, FAST exams, bladder volumes, and more. It also provides high resolution images of musculoskeletal structures and abdominal organs, such as the abdominal aorta, urinary bladder, kidneys, liver, gallbladder and abdominal free fluid. **With its small size and affordable price, Signos offers a choice that was previously not available in an ultrasound system.**



**Deliver Patient Care Quickly and Affordably with Signos Personal Ultrasound**



Signostics, Inc.  
 260 Sheridan Avenue, Suite 410  
 Palo Alto, CA 94306  
 Phone: 650-327-4000  
 Fax: 650-327-4002  
 Toll Free: 877-5-SIGNOS  
 sales@signosticsmedical.com  
 www.signosticsmedical.com